

Mapping Analyte Distributions in Surrogate Nuclear Melt Glass Using Laser-Induced Breakdown Spectroscopy and Micro X-Ray Fluorescence

LTC Michael Shattan, Ashley Stowe, Kathryn McIntosh,
John Auxier II, Christian Parigger, and Howard Hall

SCIX 2017

Reno, NV

October 2017

Outline

- Research Goals and Approach
- Motivation and Background
- Surrogate Melt Glass Production
- X-ray/Gamma Spectroscopy Comparison
- Micro-XRF Initial Results
- Laser Induced Breakdown Spectroscopy (LIBS) Initial Results
- Conclusions
- Path Forward

Research Goals:


- Explore feasibility of portable LIBS and micro-XRF systems as methods of field screening for real debris
- Develop a LIBS Capability to rapidly screen beads for production quality control
- Complete 3D elemental mapping of surrogate debris to determine uranium and other elemental migration patterns during debris formation

Scientific Approach:

- Produce melt glass “test sets” which isolate individual analytes
- Use portable COTS LIBS systems for line identification and field use
 - “Insight” benchtop LIBS by Photon Machines (Y12)
 - Nd:YAG, 266 nm, 80 mJ/pulse, 8 ns pulse width, Ocean Optics LIBS 2500 spectrometer
 - Z500-ER handheld LIBS by SciAps (Y12)
 - Nd:YAG, 1064 nm, 5 mJ/pulse, 1 ns pulse width, internal spectrometer
- Use traditional LIBS laboratory set up and physical sample sectioning for elemental mapping (UT MABE)
 - Nd:YAG, 355nm, 10ns pulse width, 80 mJ/pulse laser Princeton SpectraPro 2300i Spectrometer
- Compare results with COTS, portable Micro-XRF mapping of the same beads
 - Bruker M4 Tornado Micro-XRF (LANL)

Melt Glass Background

Nuclear Event



Type of Burst

- Environmental Analysis



Matrix Formation




Analysis Process



NTNF GCTF Tactical CONOP

- ✓ NTNF GCTF supports USG Nuclear Forensic Attribution
- ✓ Operates in same space as "Post-Det" *Consequence Management Forces*
- ✓ Key Tasks: Deploy / Route Recon / FOB Establishment / Planning / Collect Ground Debris Material



WMSC Collect



UT Surrogate Debris Formulas

- Synthetic Testable Urban Formulation for Forensics (STUFF)
 - Trinitite ; New York City; Houston
- Nuclear Under Ground Engineered Test Surrogates (NUGETS)
- Marine Urban Debris (MUD)

Soil Composition from GZ




Building Composition




City Infrastructure




Synthesis Process




Inherent Variability in Morphology and Chemical Distributions




X-Rays and Gamma Rays as Signatures for Natural Uranium



X-ray/Gamma-ray Experimental Setup










Melt Glass X-Ray/Gamma Spectrum






- Confirmation of the presence of uranium
- Long dwell time required
- No local concentration

Initial LANL Micro-XRF mapping confirms chemical heterogeneity (1 of 2)















kmb1625003

40 kV, 500 μ A
20 μ m spot
20 μ m distance
1 ms/pixel = 15 min



Initial LANL Micro-XRF mapping confirms chemical heterogeneity (2 of 2)



kmb1625204

40 kV, 500 μ A

20 μ m spot



20 μ m step

1 ms/pixel = 21 min


Melt glass production for specific LIBS line identification

Melt Glass Bead Ingredients for Line Identification

NYC Recipe	Nominal %
SiO ₂	60.48
NaO	3.23
Fe ₂ O ₃	7.63
Al ₂ O ₃	15.10
CaO	6.35
KOH	3.55
MgO	2.65
S	0.06
BaO	0.06
MnO	0.12
Ca ₃ (PO ₄) ₂	0.10
TiO	0.61
UNH	0.06
Total	100.00



Broadband LIBS Spectrum to identify interference free lines




Analytical lines for non-radioactive elements


Analytical Lines Selected

Peak (nm)	Source	A (s ⁻¹)	Eu (eV)	g
248.327	Fe I	4.80E+08	4.9913	11
279.800	Mg I	4.79E+08	8.8637	6
280.270	Mg II	2.57E+08	4.4224	2
288.158	Si I	2.17E+08	5.0823	3
328.560	Na II	1.10E+08	37.0948	5
373.486	Fe I	9.01E+07	4.1777	11
394.401	Al I	4.99E+07	4.1427	2
396.152	Al I	9.85E+07	3.1427	2
403.076	Mn I	1.70E+07	3.0751	8
440.475	Fe I	2.75E+07	4.3714	9
453.478	Ti I	6.87E+07	3.5693	9
453.557	Ti I	lot Available	3.5587	7
455.403	Ba I	1.18E+08	2.7218	4
521.262	S I	8.47E+07	17.4457	6
542.999	Fe I	6.00E+07	12.8787	10
553.548	Ba I	1.19E+08	2.2392	3
561.863	Fe I	2.24E+06	6.4149	5
585.745	Ca I	6.60E+07	5.0486	5
588.995	Na I	6.16E+07	2.1044	4
616.217	Ca I	4.77E+07	3.9104	3
656.279	H I	4.41e+07	12.0875	18
766.429	Fe I	3.37E+05	4.6076	7
766.490	K I	3.80E+07	1.6171	4
769.897	K I	3.75E+07	1.6100	2


Aluminum Analytical Lines



Barium Analytical Line



Magnesium Analytical Lines




Uranium line identification via SciAps Z500-ER Handheld LIBS

Enhanced Uranium Content Beads	
UNH Mass	Nominal %
NYC (7.2 mg)	0.06%
12.6 mg	0.11%
50.0 mg	0.42%
500.0 mg	4.17%




Uranium Line Identification


U(II) 409.013 nm




U(I) 502.738 nm



U(I) 591.539 nm




U(I) 682.691 nm



Uranium Content Variation


12.6 mg Bead

U(I) 502.738 nm




12.6 mg Bead

U(I) 591.539 nm




7.2 mg Bead

U(I) 502.738 nm




7.2 mg Bead


U(I) 591.539 nm





Experimental Setup for Chemical Mapping



Time Delay and Gate Optimization






Four Fe Lines Chosen for Analysis; One Line of Spots Ablated



Lines Used for Iron Mapping			
Source	Wavelength (nm)	Akig (s ⁻¹)	Eu (eV)
Fe (I)	422.22128	4.03E+07	5.3852074
Fe (I)	424.74253	2.13E+08	6.2864786
Fe (I)	425.83156	1.78E+04	2.9980448
Fe (I)	427.11535	1.64E+08	5.3515703


50 pulses per location




Ablation Craters

- diameter: 290-300 μm
- Depth: 300-500 μm

Positions 10-14 Provide Sample Information



Fe concentrations via peak ratios show self consistency



Conclusions

- Micro-XRF and LIBS both have the capability to map analytes in melt glass production
- Both have COTS instruments that could be used in a forward deployed environment
- The complex nature of the sample presents challenges to both techniques
- Multiple instruments may allow sample collectors greater confidence in the fidelity of their samples before shipment to offsite locations

Next Steps

- Round Robin 2D Chemical Mapping with Micro-XRF and LIBS systems on a common sample set
- 3D mapping via physical sectioning
- Increase sample throughput to measure variability in the melt glass production process

Acknowledgements

- Dr. Howard Hall (UTK, INS)
- Dr. John Auxier II (UTK, INS)
- Dr. Matthew Cook (UTK, INS)
- Dr. Graham Walford (UTK, ORNL)
- Dr. Christian Parigger (UTSI)
- Dr. Ashly Stowe (Y12)
- Cody Nizinski (UTK)
- Asia Conley (UTK)
- Jason Matheny (UTK)



 Thank You!

